


Empower® Software Tutorial: Step-by-Step Project Creation





Click on "Configure the System."



Filter By: Default

Edit View

Update

Empower 3 Configuration

	Name	Owner	Create Date	Full Audit Trail	Locked	Comments	Tablespace Quota
1	CIA_Default	System	22/02/2011 13:19:43 CET	<input type="checkbox"/>	No Lock	CIA demonstration project	
	Fields	System	22/02/2011 13:20:34 CET	<input type="checkbox"/>	No Lock	Project containing custom fields available to copy into other projects	
		System	22/02/2011 13:21:06 CET	<input type="checkbox"/>	No Lock	Default project	
		hilal	25/02/2011 12:17:20 CET	<input type="checkbox"/>	No Lock	Default project	
	FAT	S					
	FAT1	hi					
		S					
		hi					
	on_Default	S				to display bath temperature and speed as control charts	
	ult	S					
	Default	S					
	efault	System	22/02/2011 13:25:59 CET	<input type="checkbox"/>	No Lock		
			26:34 CET	<input type="checkbox"/>	No Lock		
			27:08 CET	<input type="checkbox"/>	No Lock	Demonstration Project	
			27:34 CET	<input type="checkbox"/>	No Lock	Demonstration Project	
			28:00 CET	<input type="checkbox"/>	No Lock	Demonstration Project	
			28:24 CET	<input type="checkbox"/>	No Lock	Default project for Waters Integrity system (2690-996-TMD)	
17	MS_De		28:24 CET	<input type="checkbox"/>	No Lock		
18	MVM_D		29:21 CET	<input checked="" type="checkbox"/>	No Lock	Demonstration Project	
19	MVM_D		30:59 CET	<input checked="" type="checkbox"/>	No Lock	Validation Template for Demonstration Project	
20	Pattern	System	22/02/2011 13:32:22 CET	<input type="checkbox"/>	No Lock	Chromatographic Pattern Matching	
21	PDA Default	System	22/02/2011 13:31:34 CET	<input type="checkbox"/>	No Lock	Default PDA Project	
22	PDA Default1	hilal	25/02/2011 12:18:00 CET	<input type="checkbox"/>	No Lock	Default PDA Project	
23	SQT_Alliance2489	System	22/02/2011 15:10:00 CET	<input checked="" type="checkbox"/>	No Lock	Qualification project for SystemsQT v1.10.	
24	SQT_AllianceHT2487	System	22/02/2011 15:07:00 CET	<input checked="" type="checkbox"/>	No Lock	Qualification project for SystemsQT v1.00.	
25	SysSuit_Default	System	22/02/2011 13:32:59 CET	<input type="checkbox"/>	No Lock	System suitability demonstration project	
26	SysSuit_Default1	hilal	25/02/2011 12:18:18 CET	<input type="checkbox"/>	No Lock	System suitability demonstration project	
27	ZQ_Default	System	22/02/2011 13:33:58 CET	<input type="checkbox"/>	No Lock	ZQ Demonstration Project	

Right-click with the mouse on
Projects
Then click on New and select
Project.



Filter By: Default

Edit View

Update

Empower 3 Configuration

	Name	Owner	Create Date	Full Audit Trail	Locked	Comments	Tablespace Quota
1	CIA_Default	System	22/02/2011 13:19:43 CET				
			13:34 CET			nto other projects	
			13:06 CET				
			13:20 CET				
			13:57 CET				
			13:34 CET				
			13:33 CET				
			13:46 CET				
			13:32 CET			th temperature and speed as control charts	
			13:16 CET				
			13:10 CET				
			13:59 CET				
			13:34 CET				
			13:08 CET				
			13:34 CET				
			13:00 CET				
			13:24 CET			996-TMD)	
			13:21 CET				
			13:59 CET				
			13:22 CET	<input type="checkbox"/>	No Lock	Chromatographic Pattern matching	
			13:34 CET	<input type="checkbox"/>	No Lock	Default PDA Project	
			13:00 CET	<input type="checkbox"/>	No Lock	Default PDA Project	
			13:00 CET	<input checked="" type="checkbox"/>	No Lock	Qualification project for SystemsQT v1.10.	
			13:00 CET	<input checked="" type="checkbox"/>	No Lock	Qualification project for SystemsQT v1.10.	
25	SysSuit_Default	System	22/02/2011 13:32:59 CET	<input type="checkbox"/>	No Lock	System suitability demonstration project	
26	SysSuit_Default1	hilal	25/02/2011 12:18:18 CET	<input type="checkbox"/>	No Lock	System suitability demonstration project	
27	ZQ_Default	System	22/02/2011 13:33:58 CET	<input type="checkbox"/>	No Lock	ZQ Demonstration Project	

New Project Wizard - Tablespace

Enter the amount of database tablespace to reserve for the new project.

Note: This value may be changed at any time.

TableSpace: MB

TableSpace Available: 135,19 MB

Full Audit Trail Support

Project Audit Trail Policies

Project Object	Comment	Confirm Identity
1 Method	Silent	<input type="checkbox"/>
2 Result	Silent	<input type="checkbox"/>
3 Sample	Silent	<input type="checkbox"/>
4 Deletion	Silent	<input type="checkbox"/>

Data Processing Techniques

Enable ApexTrack Integration

Default Algorithm:

< Back

Next >

Cancel

Help

You can leave everything as default and click on "Next" Or Uncheck the "Full Audit Trail" box (this option is used to record everything that happens at the project level) For the rest, leave it as default.





Filter By: Default

Edit View

Update

Empower 3 Configuration

	Name	Owner	Create Date	Full Audit Trail	Locked	Comments	Tablespace Quota
1	CIA_Default	System	22/02/2011 13:19:43 CET	<input type="checkbox"/>	No Lock	CIA demonstration project	
			13:34 CET	<input type="checkbox"/>	No Lock	Project containing custom fields available to copy into other projects	
			13:06 CET	<input type="checkbox"/>	No Lock	Default project	
			13:20 CET	<input type="checkbox"/>	No Lock	Default project	
			13:57 CET	<input checked="" type="checkbox"/>	No Lock	Full Audit Trail Demonstration Project	
			13:34 CET	<input checked="" type="checkbox"/>	No Lock	Full Audit Trail Demonstration Project	
			13:33 CET	<input type="checkbox"/>	No Lock		
						ications and to display bath temperature and speed as control charts	
			13:10 CET	<input type="checkbox"/>	No Lock		
			13:59 CET	<input type="checkbox"/>	No Lock		
			13:34 CET	<input type="checkbox"/>	No Lock		
			13:08 CET	<input type="checkbox"/>	No Lock	Demonstration Project	
			13:34 CET	<input type="checkbox"/>	No Lock	Demonstration Project	
			13:00 CET	<input type="checkbox"/>	No Lock	Demonstration Project	
			13:24 CET	<input type="checkbox"/>	No Lock	Default project for Waters Integrity system (2690-996-TMD)	
			13:21 CET	<input checked="" type="checkbox"/>	No Lock	Demonstration Project	
			13:59 CET	<input checked="" type="checkbox"/>	No Lock	Validation Template for Demonstration Project	
			13:22 CET	<input type="checkbox"/>	No Lock	Chromatographic Pattern Matching	
			13:34 CET	<input type="checkbox"/>	No Lock	Default PDA Project	
			13:00 CET	<input type="checkbox"/>	No Lock	Default PDA Project	
			13:00 CET	<input checked="" type="checkbox"/>	No Lock	Qualification project for SystemsQT v1.10.	
			13:00 CET	<input checked="" type="checkbox"/>	No Lock	Qualification project for SystemsQT v1.00.	
25	SysSuit_Default	System	22/02/2011 13:32:59 CET	<input type="checkbox"/>	No Lock	System suitability demonstration project	
26	SysSuit_Default1	hilal	25/02/2011 12:18:18 CET	<input type="checkbox"/>	No Lock	System suitability demonstration project	
27	ZQ_Default	System	22/02/2011 13:33:58 CET	<input type="checkbox"/>	No Lock	ZQ Demonstration Project	

New Project Wizard - Options



Enabled Options:

- Photo Diode Array
- System Suitability
- GPC
- GPC/V
- Mass Spectrometry
- CE/CIA
- GPC/LS
- Dissolution
- Chemical Structures

< Back

Next >

Cancel

Help

Leave by default and
click on Next.





Filter By: Default

Edit View

Update

Empower 3 Configuration

	Name	Owner	Create Date	Full Audit Trail	Locked	Comments	Tablespace Quota
1	CIA_Default	System	22/02/2011 13:19:43 CET	<input type="checkbox"/>	No Lock	CIA demonstration project	

	Name	Owner	Create Date	Full Audit Trail	Locked	Comments	Tablespace Quota
2			22/02/2011 13:34:06 CET	<input type="checkbox"/>	No Lock		
3			22/02/2011 13:57:20 CET	<input type="checkbox"/>	No Lock		
4			22/02/2011 13:57:34 CET	<input type="checkbox"/>	No Lock		
5			22/02/2011 13:57:33 CET	<input type="checkbox"/>	No Lock		
6			22/02/2011 13:46:00 CET	<input type="checkbox"/>	No Lock		
7			22/02/2011 13:32:00 CET	<input type="checkbox"/>	No Lock	use to show f1 and f2 calculations and to display bath temperature and speed as control charts	
8			22/02/2011 13:16:00 CET	<input type="checkbox"/>	No Lock	GC Demonstration Project	
9			22/02/2011 13:10:00 CET	<input type="checkbox"/>	No Lock		
10			22/02/2011 13:59:00 CET	<input type="checkbox"/>	No Lock		
11			22/02/2011 13:34:00 CET	<input type="checkbox"/>	No Lock		
12			22/02/2011 13:08:00 CET	<input type="checkbox"/>	No Lock	Demonstration Project	
13			22/02/2011 13:34:00 CET	<input type="checkbox"/>	No Lock	Demonstration Project	
14			22/02/2011 13:00:00 CET	<input type="checkbox"/>	No Lock	Demonstration Project	
15			22/02/2011 13:24:00 CET	<input type="checkbox"/>	No Lock	Default project for Waters Integrity system (2690-996-TMD)	
16			22/02/2011 13:21:00 CET	<input checked="" type="checkbox"/>	No Lock	Demonstration Project	
17			22/02/2011 13:59:00 CET	<input checked="" type="checkbox"/>	No Lock	Validation Template for Demonstration Project	
18			22/02/2011 13:22:00 CET	<input type="checkbox"/>	No Lock	Chromatographic Pattern Matching	
19			22/02/2011 13:34:00 CET	<input type="checkbox"/>	No Lock	Default PDA Project	
20			22/02/2011 13:00:00 CET	<input type="checkbox"/>	No Lock	Default PDA Project	
21			22/02/2011 13:00:00 CET	<input checked="" type="checkbox"/>	No Lock	Qualification project for SystemsQT v1.10.	
22			22/02/2011 13:00:00 CET	<input checked="" type="checkbox"/>	No Lock	Qualification project for SystemsQT v1.00.	
23			22/02/2011 13:32:59 CET	<input type="checkbox"/>	No Lock	System suitability demonstration project	
24			22/02/2011 13:32:59 CET	<input type="checkbox"/>	No Lock	System suitability demonstration project	
25	SysSuit_Default	System	22/02/2011 13:32:59 CET	<input type="checkbox"/>	No Lock	System suitability demonstration project	
26	SysSuit_Default1	hilal	25/02/2011 12:18:18 CET	<input type="checkbox"/>	No Lock	System suitability demonstration project	
27	ZQ_Default	System	22/02/2011 13:33:58 CET	<input type="checkbox"/>	No Lock	ZQ Demonstration Project	

New Project Wizard - Access Control



Allowed Access
 Owner Only
 Owner and Group
 Owner, Group and World

Select the users that should have access to this project.

Group User Type
 User's Own Type

Select the type of user access given to the group(s).

Allow Access to Groups
 Administrators
 Guests

Select the group(s) to be given access to the project.

World User Type
 User's Own Type

Select the type of access given other users.

< Back

Next >

Cancel

Help

Select Owner, Group, and World.
 Check the boxes for Administrators and Guests.
 Then click on Next.



Formation Wat...

E Pro

Defaults as Syst...



Processing Serv...

Microsoft Powe...

System/Admini...

FR



11:34



Filter By: Default

Edit View

Update

Empower 3 Configuration

	Name	Owner	Create Date	Full Audit Trail	Locked	Comments	Tablespace Quota
1	CIA_Default	System	22/02/2011 13:19:43	<input type="checkbox"/>			
2	Custom_Field	System	22/02/2011 13:20:34	<input type="checkbox"/>			
3	CET			<input type="checkbox"/>	No Lock	Full Audit Trail Demonstration Project	
33	CET			<input type="checkbox"/>	No Lock		
46	CET			<input type="checkbox"/>	No Lock		
32	CET			<input type="checkbox"/>	No Lock	use to show f1 and f2 calculations and to display bath temperature and speed as control charts	
16	CET			<input type="checkbox"/>	No Lock	GC Demonstration Project	
10	CET			<input type="checkbox"/>	No Lock		
59	CET			<input type="checkbox"/>	No Lock		
34	CET			<input type="checkbox"/>	No Lock		
08	CET			<input type="checkbox"/>	No Lock	Demonstration Project	
34	CET			<input type="checkbox"/>	No Lock	Demonstration Project	
00	CET			<input type="checkbox"/>	No Lock	Demonstration Project	
24	CET			<input type="checkbox"/>	No Lock	Default project for Waters Integrity system (2690-996-TMD)	
21	CET			<input checked="" type="checkbox"/>	No Lock	Demonstration Project	
59	CET			<input checked="" type="checkbox"/>	No Lock	Validation Template for Demonstration Project	
22	CET			<input type="checkbox"/>	No Lock	Chromatographic Pattern Matching	
34	CET			<input type="checkbox"/>	No Lock	Default PDA Project	
00	CET			<input type="checkbox"/>	No Lock	Default PDA Project	
00	CET			<input checked="" type="checkbox"/>	No Lock	Qualification project for SystemsQT v1.10.	
00	CET			<input checked="" type="checkbox"/>	No Lock	Qualification project for SystemsQT v1.00.	
59	CET			<input type="checkbox"/>	No Lock	System suitability demonstration project	
26	SysSuit_Default1	hilal	25/02/2011 12:18:18 CET	<input type="checkbox"/>	No Lock	System suitability demonstration project	
27	ZQ_Default	System	22/02/2011 13:33:58 CET	<input type="checkbox"/>	No Lock	ZQ Demonstration Project	

New Project Wizard - Copy Selection

Select items to copy:

Copy

- View Filters
- Custom Fields
- Methods
- Preferences

From Project:

- Projects
 - CIA_Default
 - Custom_Field
 - Defaults
 - Defaults1
 - Defaults_FAT
 - Defaults_FAT1

< Back

Next >

Cancel

Help

Select the Default project and click on Next. This Default project will serve as a template to create the other projects.





Filter By: Default

Edit View


Update

Empower 3 Configuration

- Projects
- Nodes
- Systems
- Libraries
- eCord
- Users
- User Groups
- User Types
- Plate Types
- System Audit Trail
- Offline System Audit Trail
- Project Archives
- Sample Archives
- Offline Project Archives
- Offline Sample Archives


E	Name	Owner	Create Date	Full Audit Trail	Locked	Comments	Tablespace Quota
1	CIA_Default	System	22/02/2011 13:00:00 CET	<input type="checkbox"/>	No Lock	CIA demonstration project	
2	Custom_Fields	System	22/02/2011 13:00:00 CET	<input type="checkbox"/>	No Lock	available to copy into other projects	
3	Defaults	System	22/02/2011 13:00:00 CET	<input type="checkbox"/>	No Lock		
4	Defaults1	hial	25/02/2011 12:00:00 CET	<input type="checkbox"/>	No Lock		
5	Defaults_FAT	System	22/02/2011 13:00:00 CET	<input type="checkbox"/>	No Lock	project	
6	Defaults_FAT1	hial	25/02/2011 12:00:00 CET	<input type="checkbox"/>	No Lock	project	
7	Demo	System	22/02/2011 13:00:00 CET	<input type="checkbox"/>	No Lock		
8	Demo1	hial	25/02/2011 12:00:00 CET	<input type="checkbox"/>	No Lock		
9	Dissolution_Default	System	22/02/2011 13:00:00 CET	<input type="checkbox"/>	No Lock	ns and to display bath temperature and speed as control charts	
10	Formation	System	21/04/2011 11:47:57 CET	<input checked="" type="checkbox"/>	No Lock	test	
11	GC_Default	System	22/02/2011 13:24:16 CET	<input type="checkbox"/>	No Lock	GC Demonstration Project	
12	GPCVLS_Default	System	22/02/2011 13:25:10 CET	<input type="checkbox"/>	No Lock		
13	GPCV_Default	System	22/02/2011 13:25:59 CET	<input type="checkbox"/>	No Lock		
14	GPC_Default	System	22/02/2011 13:26:34 CET	<input type="checkbox"/>	No Lock		
15	LIMS_ASCII_Andi	System	22/02/2011 13:27:08 CET	<input type="checkbox"/>	No Lock	Demonstration Project	
16	LIMS_LabSystems	System	22/02/2011 13:27:34 CET	<input type="checkbox"/>	No Lock	Demonstration Project	
17	LIMS_PE	System	22/02/2011 13:28:00 CET	<input type="checkbox"/>	No Lock	Demonstration Project	
18	MS_Default	System	22/02/2011 13:28:24 CET	<input type="checkbox"/>	No Lock	Default project for Waters Integrity system (2690-996-TMD)	
19	MVM_Defaults	System	22/02/2011 13:29:21 CET	<input checked="" type="checkbox"/>	No Lock	Demonstration Project	
20	MVM_Defaults_Val_Template	System	22/02/2011 13:30:59 CET	<input checked="" type="checkbox"/>	No Lock	Validation Template for Demonstration Project	
21	PatternMatch_Default	System	22/02/2011 13:32:22 CET	<input type="checkbox"/>	No Lock	Chromatographic Pattern Matching	
22	PDA Default	System	22/02/2011 13:31:34 CET	<input type="checkbox"/>	No Lock	Default PDA Project	
23	PDA Default1	hial	25/02/2011 12:18:00 CET	<input type="checkbox"/>	No Lock	Default PDA Project	
24	SQT_Alliance2489	System	22/02/2011 15:10:00 CET	<input checked="" type="checkbox"/>	No Lock	Qualification project for SystemsQT v1.10.	
25	SQT_AllianceHT2487	System	22/02/2011 15:07:00 CET	<input checked="" type="checkbox"/>	No Lock	Qualification project for SystemsQT v1.00.	
26	SysSuit_Default	System	22/02/2011 13:32:59 CET	<input type="checkbox"/>	No Lock	System suitability demonstration project	
27	SysSuit_Default1	hial	25/02/2011 12:18:18 CET	<input type="checkbox"/>	No Lock	System suitability demonstration project	
28	ZQ_Default	System	22/02/2011 13:33:58 CET	<input type="checkbox"/>	No Lock	ZQ Demonstration Project	

Congratulations, you have created your first project!




Waters


THE SCIENCE OF WHAT'S POSSIBLE.™




Guide to Creating Method Sets and Sample Sets







Select the project and the channel on which you will perform your injections. Then click OK.



2695_2487 in Formation as System/Administrator - Run Samples

File Edit View Inject Actions Customize Help



Apply Table Preferences

Sample Set Method

Run and Report

Continue on Fault

Sample Name:

Function: Inject Samples

Method Set:

Vial: 1

Develop Methods


Injection Volume: 10,0

Run Time: 10,00

Options...



Click on "Edit."



Single Samples Sample Sets R

Flow (mL/min)



Pressure (psi)



Instrument Method:

Edit


Monitor

Setup

Sample Set Time Remaining: 0:00:00:00

Total Samples Time Remaining: 0:00:00:00

New Sample Set Time: 0:00:00:00



Temperature (°C)

For Help, press F1

System Idle

2695_2487 in Formation as System/Administrator - Run Samples

File Edit View Inject Actions Customize Help



Apply Table Preferences

Sample Set Method

Run and Report

Continue on Fault

Sample Name:
Function: Inject Samples
Method Set:
Vial: 1
Injection Volume: 10.0
Run Time: 10.00



Untitled in Formation as System/Administrator - Instrument Method Editor

File Edit View Help

W2690/5 W2487

General | Degas | Events | Flow | Temperature | Solvents | Channel

General System Parameters

Stroke Volume (uL)	50uL (flow rates <= 1.23 mL/min)	<input checked="" type="checkbox"/> Bubble Detect
Syringe Draw Rate(uL/sec)	25uL (flow rates <= 0.53 mL/min)	
	50uL (flow rates <= 1.23 mL/min)	
	100uL (flow rates <= 3.030 mL/min)	
	130uL (flow rates <= 10.000 mL/min)	
Depth Of Needle(mm)	0.0	Chart Out %A
Column Position	No Change	Needle Wash Time Normal
Equilibration Time (mins)	0.00	

Ready

Introduction of the parameters for the pump instrument method

Single Samples Sample

Flow (mL/min) Pressure (psi)



Edit

Monitor

Setup

New Sample Set Time: 0:00:00:00

Sample Set

0,000 0,002 0,004 0,006 0,008 0,010
Liters

Temperature (°C)

For Help, press F1

System Idle

The screenshot displays the Waters software interface. The main window is titled "2695_2487 in Formation as System/Administrator - Run Samples". The menu bar includes "File", "Edit", "View", "Inject", "Actions", "Customize", and "Help". The toolbar contains various icons for file operations and system controls. Below the toolbar, there are input fields for "Sample Name:", "Function: Inject Samples", "Method Set:", "Vial: 1", "Injection Volume: 10,0", and "Run Time: 10,00".

An "Instrument Method Editor" window is open, titled "Untitled in Formation as System/Administrator - Instrument Method Editor". Its menu bar includes "File", "Edit", "View", and "Help". The toolbar of this window shows icons for "W2690/5" and "W2487". A red arrow points to the "W2487" icon. Below the toolbar, the "General" tab is selected, showing "Wavelength Mode" set to "Single". Under "Channel 1", the "Enable" checkbox is checked, and the "Name" is "2487Channel 1". Under "Channel 2", the "Enable" checkbox is unchecked, and the "Name" is "2487Channel 2".

At the bottom of the main window, there are buttons for "Edit", "Monitor", and "Setup". The status bar shows "System Idle" and "New Sample Set Time: 0:00:00:00". A "Sample Set" dropdown menu is visible, and a "Temperature (°C)" display is on the right.

Introduction of the parameters of the instrument method of the Detector.


2695_2487 in Formation as System/Administrator - Run Samples

File Edit View Inject Actions Customize Help

Apply Table Preferences Sample Set Method Run and Report Continue on Fault

formation MI in Formation as System/Administrator - Instrument Method Editor

File Edit View Help

- New Ctrl+N
- Open... Ctrl+O
- Save Ctrl+S
- Save As... 
- Exit

Wavelength Mode
 Single Dual

Channel 1
 Enable
Name : 2487Channel 1
Description:

Channel 2
 Enable
Name : 2487Channel 2
Description:

Save current Instrument Method

Names:
Alliance
formation MI

Name: formation MI

Method Comments:

Save Cancel Help

Save the active document with a new name

Flow (mL/min) Pressure (psi)

Single Samples Sample

System Idle

New Sample Set Time: 0:00:00:00

Sample Set 0,000 0,002 0,004 0,006 0,008 0,010 Liters



Apply Table Preferences

Sample Set Method

Run and Report

Continue on Fault

Sample Name:

Function: Inject Samples

Method Set:

Vial: 1

Develop Methods



Injection Volume: 10,0

Run Time: 10,00

Options...



Select the instrument method (MI), then click on Develop Methods

Single Samples Sample Sets

Flow (mL/min)



Pressure (psi)




Instrument Method:

formation MI


Sample Set Time Remaining: 0:00:00:00

Total Samples Time Remaining: 0:00:00:00

New Sample Set Time: 0:00:00:00




Temperature (°C)




Select the MI

Then click on "Next."



Then click on "Next."



Give a name to the Method Set (MS)
Then click on Finish.

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™



Empower® Software Essentials: Step-by-Step Tutorial for Processing Method Creation and Data Analysis

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



Mes documents Nokia PC Suite Google Earth1

Poste de travail Favoris réseau

Favoris réseau µTorrent

Corbeille Microsoft Office Word 2007

Internet Explorer Thermo.doc

Empower SOLAAR

Waters Column Selection Partie par million.docx


MassLynx V4.1 Nero StartSmart

Mozilla Firefox YRezig-CV.pdf

Empower Pro [min] [max] [close]

LOGIN LOGOUT ABOUT HELP

User: System/Administrator
DB: Local

 Run Samples	 Browse Project	 Configure System
 Process Data	 Review Data	 Print Data

Mes documents Nokia PC Suite Google Earth1

Poste de travail Favoris réseau

Favoris réseau µTorrent

Corbeille Microsoft Office Word 2007

Internet Explorer Thermo.doc

Empower SOLAAR

Waters Column Selection Partie par million.docx

MassLynx V4.1 Nero StartSmart

Mozilla Firefox YRezig-CV.pdf

Empower Pro

Browse Project

Use 'Browse Project' to open the Empower Project Window.

Select the desired project from the displayed list.

Use the Project Window to:


- Create methods for acquiring, processing and reporting.
- Copy methods and data between projects or to Windows folders on the desired drive.
- Review data.
- Preview data.
- Acquire data.

Project to browse


- Projects
 - Defaults
 - test
 - Valid_OQ1
 - Valid_OQApex

OK Cancel Help



File Edit View Plot Process Navigate Options Window Help



0 0



Name	Retention Time (min)	Area ($\mu\text{V} \cdot \text{sec}$)	% Area	Height (μV)	Int Type	Amount	Units	Peak Type	Peak Codes



Processing Method Wizard


You can use the Processing Method wizard to edit the current processing method or to create a new processing method.

Note: If you edit an existing processing method using the Processing Method wizard, the software will clear the calibration associated with the edited processing method and current 2D channel.



Create a New Processing Method
 Edit an Existing Processing Method

OK Cancel Help

Name	Retention Time (min)	Area ($\mu\text{V}^*\text{sec}$)	% Area	Height (μV)	Int Type	Amount	Units	Peak Type	Peak Codes





Select the part of the chromatogram that you want to process, either by:
1- Entering the start and end time of the chromatogram.



Integration - Peak Rejection - 486

If the integration is not correct, press the Back button to readjust the chromatogram zoom level to a different region and try again. To reject peaks based on area or height, select the smallest peak of interest and click the Minimum Area or Minimum Height box.




Minimum Area 10000 Minimum Height 0

Test


< Précédent Suivant > Annuler

Name	Retention Time (min)	Area (μV*sec)	% Area	Height (μV)



File Edit View Plot Process Navigate Options Window Help

37,50 37,50 2,000 2,000 100000 0



Processing Method Name - 486

Method Name:

Click Finish to save the modified processing method and return to Review where the saved processing method will be used to process the current chromatogram.


< Précédent **Terminer** Annuler

Name	Retention Time (min)	Area (µV*sec)	%
1 Acetone	1,249	198704	2
2 Acetophenone	2,132	196045	27,71 39424 BB 2,500 Found
3 Propiophenone	3,110	155486	21,98 26228 BB 2,500 Found
4 Butyrophenone	4,781	157316	22,23 20122 BB 2,500 Found

Enter a name example:
Training

File Edit View Plot Process Navigate Options Window Help

37,50 37,50 2,000 2,000 100000 0



Name	Retention Time (min)	Area (µV*sec)	% Area	Height (µV)	Int Type	Amount	Units	Peak Type	Peak Codes
Acetone	1,249	198704	28,08	47020	BV	937,500		Found	
Acetophenone	2,132	196045	27,71	39424	BB	2,500		Found	
Propiophenone	3,110	155486	21,98	26228	BB	2,500		Found	
Butyrophenone	4,781	157316	22,23	20122	BB	2,500		Found	

Vial	Inj Vol (uL)	# of Injs	Label	SampleName	Sample Type	Level	Method Set / Report Method	SampleWeight	Dilution	Altered
1	20.0	1	S0101	PQ Std. 2.5x	Standard	1	LC Demo Method Set	1,00000	1,00000	<input checked="" type="checkbox"/>
2	20.0	1	S0102	PQ Std 5.0x	Standard	2	LC Demo Method Set	1,00000	1,00000	<input checked="" type="checkbox"/>
3	20.0	1	S0103	PQ Std 10x	Standard	3	LC Demo Method Set	1,00000	1,00000	<input checked="" type="checkbox"/>
4	20.0	1	U0101	PQ Unk. 1	Unknown		LC Demo Method Set	1,00000	1,00000	<input checked="" type="checkbox"/>
5	20.0	1	U0102	PQ Unk. 2	Unknown		LC Demo Method Set	1,00000	1,00000	<input checked="" type="checkbox"/>
6	20.0	1	U0103	PQ Unk. 3	Unknown		LC Demo Method Set	1,00000	1,00000	<input checked="" type="checkbox"/>
7	20.0	1	U0104	PQ Unk. 4	Unknown		LC Demo Method Set	1,00000	1,00000	<input checked="" type="checkbox"/>

Component Editor

File Edit View Help

SampleSet Type: STANDARDS ONLY

Current Vial
Row: 2 Vial: 2 Level: 2 Type: Standard


Components

Component	Value	Units (Vial)


Current All Samples /

Prev Next OK Cancel

"Click on Edit."





Click on "Copy Component From Process Method."



Select the previously created Process Method and click on Open







To perform the process for unknowns and determine their concentration:

1. First, select the Standards in ascending order of concentration (from the lowest to the highest).
2. Then, choose the unknown samples.
3. Finally, click on the "Process" icon.

Defaults as System/Administrator - Project

File Edit View Tools Database Help

Filter By: Default

SampleName	Vial	Injection	Sam
1 PQ Unk. 4	7	1 Unkn	
2 PQ Unk. 3	6	1 Unkn	
3 PQ Unk. 2	5	1 Unkn	
4 PQ Unk. 1	4	1 Unkn	
5 PQ Std 10x	3	1 Stand	
6 PQ Std 5.0x	2	1 Stand	
7 PQ Std. 2.5x	1	1 Stand	

Defaults - Background Processing and Reporting

Processing

Process

Use acquisition method set (i.e. from the sample set used to acquire data)

Use specified method set

Use specified processing method

Clear Calibration

LC PQ
LC PQ ApexTrack
pur2
TEST
test MP

Reporting

Print Fax

Use acquisition method set (i.e. from the sample set used to acquire data)

Use specified method set

Use specified report method

Exporting

Export

Use acquisition method set (i.e. from the sample set used to acquire data)

Use specified method set


Use specified export method

OK Cancel Help

7 total


Formation Waters Microsoft PowerPoi... E Pro Defaults as System/... Calculator FR 10:08

- 1- Click on "Use Specified Processing Method."
- 2- Select the processing method that you have created.
- 3- Click on "Clear Calibration."
- 4- Click on "ok."



Name	Retention Time (min)	Area (μV*sec)	% Area	Height (μV)	Int Type	Amount	Units	Peak Type	Peak Codes
1 Acetone	1,249	802579	28,00	187798	BV	402,970		Found	Q09
2 Acetophenone	2,133	794688	27,73	157462	VV	10,099		Found	Q09
3 Propiophenone	3,111	631403	22,03	105414	VV	10,132		Found	Q09
4 Butyrophenone	4,781	637338	22,24	81063	VB	10,064		Found	Q09

Here are the concentrations of each peak.



By clicking on the "Calibration Curve" icon:

- 1- You obtain the calibration curve generated by the three standards.
- 2- The equation of the straight line.
- 3- As well as the R² value.

Defaults as System/Administrator - Project

File Edit View Tools Database Application Help

Filter By: Default

Sample Sets

SampleName Vi

1 PQ Unk. 1 4

2 PQ Unk. 2 5

3 PQ Unk. 3 6

4 PQ Unk. 4 7

5 PQ Std. 2.5x 1

6 PQ Std 5.0x 2

7 PQ Std 10x 3

8 PQ Std 5.0x 2

9 PQ Unk. 1 4

10 PQ Unk. 3 6

11 PQ Unk. 4 7

12 PQ Std. 2.5x 1

13 PQ Unk. 2 5

14 PQ Std 10x 3

Chromatograms

Zeroed Baseline Chromatograms

Tables

LC Calibration Summary

LC Calibration Summary

LC Calibration Summary

Instrument Methods

Processing Methods

Report Methods

Export Methods

Method Sets

Sample Set Methods

System Information

Sample Sets

Sample Components

Component Status

Drawing Objects

Special Information

Fields

Result Custom Views

Composite Groups

Page Breaks

Acquisition Logs

Internal Standard Groups

Open Report Method

Please select the Report Method that you would like to use to preview the data that you have selected:

Use the Report Method Default Individual Report in the acquisition Method Set LC Demo Method Set.

Use the Report Method named Default.

Use a Report Method that was generated to be appropriate for the selected data.

Use the following Report Method: Chromatogram Labels Report

Use the currently open Report Method.

Chromatogram Labels Report

Component Summary

Default Individual Report

Injection Summary Report

LC Calibration Report

Multi Sample Summary

Overlay Report

Peak Summary Report

Replicate Injection Report

OK

For Help, press F1

Project: Defaults

Individual

Creating a report:
1. Click on "Use the following Report Method."

2. Select one of the reports from the list.

3. Click on "OK."
In this example, we will use the "Default Individual Report."




1

2



Default Individual Report

SAMPLE INFORMATION			
Sample Name:	PQ Vial 1	Acqfile by:	System
Sample Type:	Unknown	Sample Set Name:	PQ Sample Set
Vial:	4	Acq. Method Set:	LC Demo Method Set
Injection #:	1	Process Inq Method:	Formation MP
Injection Volume:	20.00 µl	Channel Name:	486
Run Time:	6.0 Minutes	Proc. Chnl. Descr.:	254µm
Date Acquired:	17/09/1997 17:17:07 EDT		
Date Processed:	02/09/2009 09:56:59 CET		



Peak Name	RT	Area	% Area	Height	Amount
1 Acetone	1.249	802579	28.00	187796	402,970
2 Acetone	2.133	794688	27.73	167462	10,099
3 Propylene	3.111	631403	22.03	106414	10,132
4 Ethylene	4.781	637338	22.24	81063	10,064


Reported by User: System
Report Method: Default Individual Report
Report Method ID: 2222
Page: 1 of 1

Project Name: Defaults
Date Printed: 02/09/2009
10:10:11 AM/ka/Angle

- 1- To print the report, click on the "Print" icon.
- 2- To edit the report and add other parameters, click on "Quit."

- Chromatograms
- Zeroed Baseline Chromatograms
- Tables
 - All Peaks Table
 - Amount Component Summary
 - Area Component Summary
 - Basic LC Peaks Table
 - Component Table
 - EP Sys Suit Result Table
 - EP Sys Suit Summary Table
 - GPC Average System Sigma
 - GPC Results Table
 - GPC Summary Results Table
 - GPC System Sigma
 - GPCV Results Table
 - Group Peaks Table
 - JP Sys Suit Result Table
 - JP Sys Suit Summary Table
 - MS Results Table
 - Pattern Match Interval Table
 - PDA Results Table
 - Peak Labels Table
 - Peak Results Table
 - Sample Set Table
 - Sign Off Table
 - Summary Table
 - Unknown Peaks Table
 - USP Sys Suit Result Table
 - USP Sys Suit Summary Table
- LC Calibration Curves
- LC Calibration Point Tables
- LC Calibration Avg. Point Tables
- Instrument Methods
- Processing Methods
- Report Methods
- Export Methods
- Method Sets
- Sample Set Methods
- System Information
- Sample Sets

SAMPLE INFORMATION			
Sample Name:	PO Unk.1	Acquired By:	Sackm
Sample Type:	Unknown	Sample Set Name:	PO Sample Set
Inst:	1	Acq. Method Set:	LC Demo Method Set
Injection #:	1	Processing Method:	Formation MP
Injection Volume:	20.00 ul	Channel Name:	486
Run Time:	6.0 Minutes	Proc. Chnl. Descr.:	25um
Date Acquired:	17.09.1997 17:17:07 EDT		
Date Processed:	02.09.2009 09:56:59 CET		




Peak Name	RT	Area	% Area	Height	Amount
1 Acetone	1,249	802579	28,00	187798	402,97C
2 Acetone	2,133	794688	27,73	157462	10,09E
3 Propylene	3,111	631403	22,03	105414	10,13C
4 Butylene	4,781	637338	22,24	81063	10,06A

We get this window.

- Chromatograms
- Zeroed Baseline Chromatograms
- Tables
 - All Peaks Table
 - Amount Component Summary
 - Area Component Summary
 - Basic LC Peaks Table
 - Component Table
 - EP Sys Suit Result Table
 - EP Sys Suit Summary Table
 - GPC Average System Sigma
 - GPC Results Table
 - GPC Summary Results Table
 - GPC System Sigma
 - GPCV Results Table
 - Group Peaks Table
 - JP Sys Suit Result Table
 - JP Sys Suit Summary Table
 - MS Results Table
 - Pattern Match Interval Table
 - PDA Results Table
 - Peak Labels Table
 - Peak Results Table
 - Sample Set Table
 - Sign Off Table
 - Summary Table
 - Unknown Peaks Table
 - USP Sys Suit Result Table
 - USP Sys Suit Summary Table
- LC Calibration Curves
- LC Calibration Point Tables
- LC Calibration Avg. Point Tables
- Instrument Methods
- Processing Methods
- Report Methods
- Export Methods
- Method Sets
- Sample Set Methods
- System Information
- Sample Sets

SAMPLE INFORMATION			
Sample Name:	PQ 1146_1	Acquired By:	Sachin
Sample Type:	Unknown	Sample Set Name:	PQ Sample Set
Unit:	1	Assay Method Set:	LC Demo Method Set
Injection #:	1	Processing Method:	Format: MP
Injection Volume:	20.00 ul	Channel Name:	695
Run Time:	6.0 Minutes	Proc. Chnl. Descr:	CS44m
Date Acquired:	17/09/1997 17:17:07 EDT		
Date Processed:	06/09/2009 09:56:59 CET		




#	Peak Name	RT	Area	% Area	Height	Amount
1	Acetone	1.249	802579	28.00	187198	402.97C
2	Acetophenone	2.133	794688	27.73	157462	10.096
3	Propiophenone	3.111	631403	22.03	105414	10.13C
4	Butyrophenone	4.781	637338	22.24	81063	10.064

System Suitability Separation Results						
#	Name	RT	USP Resolution	USP Tailors	K Prime	Width @ Tancent
1	Acetophenone	2.133				
2	Propiophenone	3.111				
3	Butyrophenone	4.781				

To add one or more new tables, simply double-click on one of the tables from the list located to the right of this window.

- Chromatograms
- Zeroed Baseline Chromatograms
- Tables
 - All Peaks Table
 - Amount Component Summary
 - Area Component Summary
 - Basic LC Peaks Table
 - Component Table
 - EP Sys Suit Result Table
 - EP Sys Suit Summary Table
 - GPC Average System Sigma
 - GPC Results Table
 - GPC Summary Results Table
 - GPC System Sigma
 - GPCV Results Table
 - Group Peaks Table
 - JP Sys Suit Result Table
 - JP Sys Suit Summary Table
 - MS Results Table
 - Pattern Match Interval Table
 - PDA Results Table
 - Peak Labels Table
 - Peak Results Table
 - Sample Set Table
 - Sign Off Table
 - Summary Table
 - Unknown Peaks Table
 - USP Sys Suit Result Table
 - USP Sys Suit Summary Table
- LC Calibration Curves
- LC Calibration Point Tables
- LC Calibration Avg. Point Tables
- Instrument Methods
- Processing Methods
- Report Methods
- Export Methods
- Method Sets
- Sample Set Methods
- System Information
- Sample Sets

SAMPLE INFORMATION			
Sample Name:	PQ Unk 1	Acquired By:	System
Sample Type:	Unknown	Sample Set Name:	PQ Sample Set
Unit:	1	Acq. Method Set:	LC Demo Method Set
Injection #:	1	Processing Method:	Formation MP
Injection Volume:	20.00 ul	Channel Name:	496
Run Time:	6.0 Minutes	Proc. Chnl. Descr:	CS1km
Date Acquired:	17/09/1997 17:17:07 EDT		
Date Processed:	0006/0009/09:56:59 CET		




Peak Name	RT	Area	% Area	Height	Amount
1 Acetone	1,249	802579	28,00	187198	402,970
2 Acetone	2,133	794688	27,73	157462	10,099
3 Propylene	3,111	631403	22,03	105414	10,132

System Suitability Generation Results				
Name	RT	USP Resolution	USP Tailing	Width at Target
1 Cut				
2 Copy				
3 New Table...				
4 Open...				

- Print...
- Print Preview
- Move To Back Ctrl+M
- Table Properties... Alt+Enter

You can also modify the tables to remove or add a parameter or a criterion:

Right-click on the table and then click on "Table Properties."



Select the parameters of your choice, then click to add to the table. Finally, click OK.

File Edit View Layout Help

100%

Chromatograms

Zeroed Baseline Chromatograms

Tables

- All Peaks Table
- Amount Component Summary
- Area Component Summary
- Basic LC Peaks Table
- Component Table
- EP Sys Suit Result Table
- EP Sys Suit Summary Table
- GPC Average System Sigma
- GPC Results Table
- GPC Summary Results Table
- GPC System Sigma
- GPCV Results Table
- Group Peaks Table
- JP Sys Suit Result Table
- JP Sys Suit Summary Table
- MS Results Table
- Pattern Match Interval Table
- PDA Results Table
- Peak Labels Table
- Peak Results Table
- Sample Set Table
- Sign Off Table
- Summary Table
- Unknown Peaks Table
- USP Sys Suit Result Table
- USP Sys Suit Summary Table

LC Calibration Curves

LC Calibration Point Tables

LC Calibration Avg. Point Tables

Instrument Methods

Processing Methods

Report Methods

Export Methods

Method Sets

Sample Set Methods

System Information

Sample Sets

SAMPLE INFORMATION

Sample Name:	PO Unk.1	Bottled By:	System
Sample Type:	Unknown	Sample Set Name:	PO Sample Set
Lot:	1	Acq. Method Set:	LC Demo Method Set
Injection #:	1	Processing Method:	Format01.MP
Injection Volume:	20.00 ul	Channel Name:	186
Run Time:	6.0 Minutes	Proc. Ch11 Descr:	254um
Date Acquired:	17/09/1997 17:17:07 EDT		
Date Processed:	10/09/2009 09:56:59 CET		

Peak Name	RT	Area	% Area	Height	Amount
1 Acetone	1,249	802579	29,00	187198	402,970
2 Acetone	2,133	794688	27,73	157462	10,099
3 Propylene	3,111	631403	22,03	105414	10,132

System Suitability Separation Results


Name	RT	USP Resolution	USP Tailing	K Prime	Width @ Tangent	% Height
1 Acetophenone	2,133					29.01
2 Propylene	3,111					19.02
3 Butylene	4,781					15.24
4 Acetone	1,249					29.32

We obtain the new table with the new parameters added. Click on 'Print Preview'



Sample Name: PQ Unk. 1 Report By: System
Sample Type: Unknown Sample Set Name: PQ Sample Set
Vial: 4 Acq. Method Set: LC Demo Method Set
Injection #: 1 Processing Method: Formation MP
Injection Volume: 20.00 ul Channel Name: 486
Run Time: 6.0 Minutes Proc. Chnl. Descr: 254nm

Date Acquired: 17/09/1997 17:17:07 EDT
Date Processed: 02/09/2009 09:56:59 CET



Peak Name	RT	Area	% Area	Height	Amount
1 Acetone	1,249	802579	28,00	187798	402,970
2 Acetophenone	2,133	794688	27,73	157462	10,099
3 Propiophenone	3,111	631403	22,03	105414	10,132
4 Butyophenone	4,781	637338	22,24	81063	10,064

Reported by: User: System Project Name: Defaults
Report Method: Default Individual Report Date Printed: 02/09/2009 09:56:59 CET

To view the rest of the report, click on "next page."

Thank

you